Answer all the questions (25 X 1 = 25)

Reference electrode and working electrode used in the estimation of Fe(II) 1 ions by potentiometry are*	point
Glass electrode and Platinum electrode	
O Platinum electrode and Standard Calomel Electrode	
Standard Calomel Electrode and Glass electrode	
Standard Calomel Electrode and Platinum electrode	

Viscosity is due to*	1 point
Resistance to fluid motion	
O Roughness of the surface	
The pressure difference between the two fluids	
O Potential energy stored in fluid	

In Mohr's method the solution needs to be near neutral, because * 1 point
Silver precipitates at low pH
Chromate forms H2CrO4 at low pH, which delays the formation of the precipitate.
Silver chloride forms at high pH,
O Potassium chromate dissolves at high pH.
The pH of a liquid solution is a measure of * 1 point
Hydrogen ion activity
O Dissolved salt content
Hydroxyl ion molarity
C Electrical conductivity
The electrolyte solution within the glass electrode (reference) of the pH 1 point meter is *
O Dilute NaCl
O Dilute HCI
Concentrated HCI
Saturated KCI

What is the working principle of conductometry? *	1 point
Measurements of pH	
Measurement of emf.	
Measurement of potential.	
Measurement of conductivity of solution.	
When sodium hydroxide is added to HCI, the H+ ions are replaced by*	1 point
Slow moving Na+ ions	
Slow moving OH- ions	
Fast moving OH- ions	
Fast moving Na+ ions	
When a strong base is added to a strong acid after the neutralization point *	1 point
Conductance remains constant	
Onductance decreases	
Conductance decreases initially and then increases gradually	
Conductance increases	

Temporary hardness of water is caused due to the presence ofdissolved solids *	1 point
Calcium bicarbonates only	
Magnesium bicarbonates only	
Sulphates and Chlorides of Calcium or Magnesium	
Calcium hydrogen carbonates and Magnesium hydrogen carbonates	
Conductivity of a solution is directly proportional to *	1 point
Current density	
Volume of the solution	
Number of ions	
Dilution	
Estimation of Fe(II) ions by potentiometry istitration. *	1 point
Complexometric	
O Acid-base	
Precipitation	
Redox	

Hardness of water is usually expressed in terms of equivalent 1 point amount of hardness. *
O K2C03
○ CaCO3
MgC03
Na2C03
Which type of reaction occurs in the following reaction AgNO3 + NaCl \rightarrow 1 point AgCl + NaNO3? *
O placement reaction
Double displacement reaction
Single replacement
Decomposition
A plot of hsp / C (reduced viscosity) vs C is a for dilute polymer 1 point solutions *
"V" shape curve
"S" shape curve
Straight line
○ Triangle

	If the ion size decreases in solutions then *	1 point
	O Does not affect the conductance	
	Conductance increases	
	Conductance decreases	
	First decreases and then increases	
	The equivalent weight of Sodium Carbonate [Na2CO3] is *	1 point
	O 40	
	55.85	
	53	
	O 63	
	Phenolphthalein color in basic medium is * Orange Steel Blue Yellow Pink	1 point
	Estimation of chloride reaction is a *	1 point
	Catalytic reaction	
	Equlibrium reaction	
	Redox reaction	
!	Precipitation reaction	

Indicator used for the estimation of hardness in EDTA method is*	1 point
Eriochrome Black – T	
O Phenolphthalein	
O Potassium dichromate	
Methyl orange	
A buffer solution comprises of *	1 point
a weak acid in solution	
a weak base in solution	
a weak acid and its conjugate base in solution	
a strong acid in solution	
In the Estimation of Fe(II) ions by potentiometry, K2Cr2O7 acts as *	1 point
Reducing agent	
Catalyst	
Oxidizing agent	
Indicator	

Volume of different concentrations (0.1, 0.2, 0.3, 0.4 and 0.5 %) of polymer 1 solution used for each viscosity measurement are *	point
varies with respect to polymer used	
remains fixed	
varies with respect to the size of the Ostwald viscometer	
varies with respect to concentration	
The unit of Conductance is *	point
O Ampere	
O Volts	
Mho	
Ohm	
Which one of the following is NOT a unit for hardness? *	point
O Degree French	
O Degree clarke	
Degree centigrade	
O Parts per million	

When a mixture of sodium carbonate and sodium hydroxide solution is titrated against HCl solution, the Phenolphthalein end point correspond to *	nt
Neutralization of hydroxyl ions and carbonate ions	
Neutralization of hydroxyl ions only	
Neutralization of carbonate ions only	
Neutralization of hydroxyl ions and half of carbonate ions	

Send me a copy of my responses.

Back

Submit

Never submit passwords through Google Forms.

This form was created inside of SRM Institute of Science and Technology. Report Abuse

Google Forms